

**Experiment 5**  
**The Determination of Acid Content in Vinegar**

Reading assignment: Chang, Chemistry 9<sup>th</sup> edition, pp. 127-131, 142-148.

**Goals**

We will use a titration to determine the concentration of acetic acid in a sample of vinegar in order to become familiar with acid-base reactions.

**Equipment and Supplies**

Digital analytical balance, 100 mL beakers (2), 10 mL graduated cylinder, 100 mL volumetric flask, 10 mL graduated pipette, pipette pump, 125 mL Erlenmeyer flasks (3), 50 mL buret, phenolphthalein indicator, ~0.1M sodium hydroxide, vinegar.


Safety Note: Safety glasses are required when performing this experiment.

**Discussion**

Acetic acid (HC<sub>2</sub>H<sub>3</sub>O<sub>2</sub>) is the active ingredient in vinegar and is responsible for its sour taste. Acetic acid is an example of a weak acid. For a 0.1 mol/L solution of acetic acid only about 1% of the acid ionizes. Compare this to a strong acid like hydrochloric acid. Very close to 100% of hydrochloric acid ionizes. A few examples of strong and weak acids are shown below:


strong acids			weak acids		
perchloric	HClO <sub>4</sub>	chloric HClO <sub>3</sub>	chlorous	HClO <sub>2</sub>	hypochlorous HClO
sulfuric	H <sub>2</sub> SO <sub>4</sub>		hydrocyanic	HCN	
hydrochloric	HCl		sulfurous	H <sub>2</sub> SO <sub>3</sub>	
hydrobromic	HBr		hydrosulfuric	H <sub>2</sub> S	
hydroiodic	HI		hydrofluoric	HF	
nitric	HNO <sub>3</sub>		nitrous	HNO <sub>2</sub>	
			phosphoric	H <sub>3</sub> PO <sub>4</sub>	
			<b>acetic</b>	<b>HC<sub>2</sub>H<sub>3</sub>O<sub>2</sub></b>	

Only one of the hydrogen atoms of the acetic acid molecule is acidic:




The structural formula for acetic acid is shown to the right.

The hydrogen attached to the oxygen atom is acidic while the other hydrogen atoms are not.




Another way of representing the acidity of acetic acid is to show its reaction with water:



Here the acetic acid protonates (transfers a proton to) the water molecule. In fact, the hydrogen ion (H<sup>+</sup>) is very reactive and doesn't exist in water. However, there is evidence that the hydronium ion does exist. So the second equation is more accurate. Sometimes equation 1 is used because of its simplicity.

To determine the amount of acetic acid in vinegar (typically 4-5% by mass) we will use an acid-base titration (neutralization reaction). In this experiment we titrate acetic acid with sodium hydroxide (a strong base). The reaction of acetic acid with sodium hydroxide is shown below:



In the reaction between acetic acid and sodium hydroxide, the acetic acid donates a proton to the hydroxide ion and is considered an acid. The hydroxide ion accepts a proton and is considered a base. The stoichiometric relationship between acetic acid and sodium hydroxide is 1:1 (from equation 2), so if the moles of NaOH used to titrate acetic acid are known, then the moles of acetic acid in a sample can be determined.


The reaction of acetic acid with sodium hydroxide is impossible to see, so we will add an indicator to determine the equivalence point for the titration.

A common acid-base indicator is phenolphthalein.

Phenolphthalein is colorless in acidic solution, but changes to pink when the solution becomes basic. At the equivalence point, the moles of acetic acid are equal to the moles of sodium hydroxide, and the solution will turn pink as more sodium hydroxide is added. Phenolphthalein is itself an acid. So when enough sodium hydroxide has been added to the acetic acid, such that there is no longer any unreacted acetic acid, the sodium hydroxide will then begin reacting with phenolphthalein. The reaction of phenolphthalein with sodium hydroxide results in a pink solution.

Buret containing sodium hydroxide solution

Erlenmeyer flask with vinegar and indicator



The picture shows a titration apparatus. The experiment makes use of a buret that contains the sodium hydroxide solution. An Erlenmeyer flask contains the vinegar solution and a couple of drops of the indicator (phenolphthalein). The sodium hydroxide solution is poured into the vinegar solution one drop at a time. The volume of sodium hydroxide required to react with all of the acetic acid in the vinegar is measured from the buret.

### Sample Calculation

A 25.00 mL sample of a hydrofluoric acid, a monoprotic acid, is titrated with a 0.155 M solution of sodium hydroxide. It is found that the indicator changes color once 31.25 mL of the sodium hydroxide solution has been added to the acid. What is the molarity of the acid?

The equation for the reaction is  $\text{HF}(\text{aq}) + \text{NaOH}(\text{aq}) \rightarrow \text{NaF}(\text{aq}) + \text{H}_2\text{O}(\text{l})$

At the equivalence point the moles of the acid (HF) are equal to the moles of the base (NaOH).\* We can use the known concentration and measured volume of the sodium hydroxide to find the number of moles of hydroxide used in the titration:

moles NaOH =  $\left( \frac{0.155 \text{ mol NaOH}}{1 \text{ L}} \right) (0.03125 \text{ L NaOH}) = 0.00484 \text{ mol NaOH}$  where the volume has been converted to liters.

At the equivalence point: mol HF = mol NaOH so moles of HF = 0.00484 mol

The concentration of the acid is equal to the number of moles divided by its volume:

$$\text{molarity of acid} = \frac{0.00484 \text{ mol HF}}{0.02500 \text{ L}} = \mathbf{0.194 \text{ mol/L}}$$

\*The solution changes color once all of the acid has reacted with the base and some of the indicator reacts with the base. So the change in color often occurs at what is called the endpoint, a little beyond the equivalence point.

**SAFETY PRECAUTIONS**

Acetic acid and sodium hydroxide are both irritants. Safety glasses must be worn at all times during this experiment. Students perform titrations individually.

1. Obtain about 60 mL of sodium hydroxide (NaOH) in a clean dry beaker. Be sure to write down the concentration of the sodium hydroxide from the bottle (or from the board at the front of the room).
2. Prepare the buret by rinsing the buret two to three times with distilled or deionized water and then with about 5 mL of sodium hydroxide.
3. Fill the buret with sodium hydroxide. You do not have to fill the buret to exactly 0.00 mL. Anywhere between 0 and 1 mL will work. Open the valve on the buret and allow a couple milliliters to flow into a waste beaker so that any air bubbles can be forced through. Record the volume of the sodium hydroxide in the buret on the data sheet.
4. The vinegar we will be using has an acetic acid concentration that would require a large volume of sodium hydroxide for this titration. So we will dilute the vinegar with water before performing the titration. Obtain 15-20 mL of vinegar in a clean and dry 100 mL beaker. Record the brand name of the vinegar.
5. Using a 10 mL volumetric pipette carefully pipet 10 mL of vinegar into a 100 mL volumetric flask.
6. Fill the flask to the 100 mL mark using distilled or deionized water. Use this diluted vinegar to perform the titrations.
7. Measure the mass of 3 clean 125 mL Erlenmeyer flasks (individually) on the digital analytical balance to the nearest 0.1 mg (0.0001 g). Record each mass on the data sheet.
8. Using a 10 mL graduated cylinder add 10 mL of the diluted vinegar into each of the weighed Erlenmeyer flasks. Record these mass on the data sheet.
9. Measure the mass of each of the Erlenmeyer flasks containing the diluted vinegar. Record these masses on the data sheet.
10. Add two drops of phenolphthalein indicator to the Erlenmeyer flask and begin the titration by slowly adding the sodium hydroxide solution from the buret to the flask. The solution should initially be clear, with no pink tint.
11. Perform the titration with sodium hydroxide and diluted vinegar. Be sure to swirl the solution in the flask while adding the sodium hydroxide solution. The endpoint for this titration is a very faint pink color that persists for more than 15 seconds. Darker pink colors indicate that you have added too much sodium hydroxide. Read the volume to 0.02 mL.
12. Perform the titration a minimum of three times until you have three sodium hydroxide volumes that are within 0.20 mL of each other. You can perform additional titrations if your delivered volumes are not close to one another.
13. Waste can be disposed of in the sink. Clean your work area before beginning calculations.

---

## Calculations

Perform these calculations on a separate sheet of paper.

1. Determine the number of moles of sodium hydroxide required to titrate the vinegar for each titration from the known molarity and the titration volume ( $V = V_2 - V_1$ ) of sodium hydroxide. Be sure that the volume of the sodium hydroxide has been converted from milliliters to liters.

$$\text{moles NaOH} = (V_{\text{NaOH}})(M_{\text{NaOH}})$$

2. The moles of acetic acid are equal to the moles of sodium hydroxide at the equivalence point.

moles acetic acid = moles sodium hydroxide

3. Determine the mass of acetic acid present in each titration from the molar mass (sometimes called molecular weight) of acetic acid and moles of acetic acid.

$$\text{molar mass} = \frac{\text{mass}}{\text{moles}} \quad \text{mass acetic acid} = (\text{molar mass})(\text{moles})$$

4. The mass of the vinegar for each titration is found from the measurements using the analytical balance ( $m_2 - m_1$ ).

5. The percentage mass of acetic acid in the vinegar is found from the mass of acetic acid and the mass of vinegar.

$$\% \text{mass} = \frac{\text{mass of acetic acid}}{\text{mass of vinegar}} \times 100$$

6. Calculate the average percent mass from your two closest values of mass percentage.

7. The vinegar was diluted by a factor of 10 in this experiment. So the mass percentage of the diluted vinegar should be multiplied by ten to find the mass percentage of the undiluted vinegar.

---

## Questions

Answer these questions on a separate piece of paper and include them with your report.

1. Show by calculation that taking 10 mL of 0.08 mol/L of acetic acid solution and diluting with water to 100 mL total volume gives a concentration that is 0.008 mol/L acetic acid.

2. Explain in three sentences why phenolphthalein changes color during the titration.

3. The color change at the endpoint should persist for at least fifteen seconds. Why is this important?

4. Describe the difference between endpoint and equivalence point.

5. Which solution would you expect to form more  $\text{H}^+$  ions in water, 0.5 M acetic acid or 0.1 M hydrochloric acid? Explain.

6. In the procedure you used, the mass of acetic acid solution was measured using a balance. Suppose we didn't have access to a balance but knew the density of the acetic acid solution. Explain how you could determine the mass of the acetic acid solution used in a titration.

---

**Observations and Notes**  
**Experiment 5**  
**The Determination of Acid Content in Vinegar**

Date \_\_\_\_\_

---

**Data Sheet**  
**Experiment 5**  
**Determination of Acid Content in Vinegar**

Name \_\_\_\_\_

Concentration of NaOH \_\_\_\_\_ mol/L

Brand name of vinegar \_\_\_\_\_

	<b>trial 1</b>	<b>trial 2</b>	<b>trial 3</b>
mass of empty Erlenmeyer flask ( $m_1$ )	_____	_____	_____
mass of flask and vinegar ( $m_2$ )	_____	_____	_____
mass of vinegar	_____	_____	_____
initial buret volume ( $V_1$ )	_____	_____	_____
final buret volume ( $V_2$ )	_____	_____	_____
volume of NaOH delivered	_____	_____	_____
moles of NaOH used	_____	_____	_____
moles of acetic acid in sample	_____	_____	_____
mass of acetic acid	_____	_____	_____
mass percent	_____	_____	_____
average mass percent	_____	_____	_____

