EXPERIMENT 8 The Preparation and Properties of Oxygen

Objective

We will observe the thermal decomposition of several oxygen-containing compounds called *oxides*. We will collect molecular oxygen and demonstrate its reactivity with several chemical elements.

* * * SAFETY GLASSES ARE REQUIRED FOR THIS EXPERIMENT. * * *

Equipment and Materials

Solid samples of sodium nitrate, sand (silicon dioxide), lead (IV) oxide, magnesium oxide, and potassium chlorate in test tubes. Potassium chlorate mixed with manganese dioxide in a large test tube, oxygen collection apparatus with three collection jars, solid phosphorus, solid sulfur, magnesium ribbon, universal indicator, and phenolphthalein.

Discussion

Oxygen is a plentiful and highly reactive element on the earth. Oxygen comprises about 21% of the earth's atmosphere. In its elemental form oxygen is a gas and *diatomic*.

Major	Components in	า Earth's	Atmosphere

Gas	Percent Mass	Gas	Percent Mass
nitrogen	78.084	neon	0.00182
oxygen	20.947	helium	0.00052
argon	0.934	krypton	0.00010
carbon dioxide	0.033		

Earth's crust contains a large amount of oxygen as well, primarily in the form of *silicates*.

Elemental Abundances in Earth's Crust

Element	Percent Mass	Element	Percent Mass
oxygen	46.6	sodium	2.8
silicon	27.7	potassium	2.6
aluminum	8.1	magnesium	2.1
iron	5.0	others	1.5
calcium	3.6		

Many <u>metals</u> react with oxygen to form compounds called **oxides**. For example, iron reacts with oxygen to form iron (III) oxide (rust): $4Fe(s) + 3O_2(g) \rightarrow 2Fe_2O_3(s)$

Metalloidswill also react with oxygen. Silicon forms silicon dioxide, the main componentof sand, when it reacts with oxygen: $Si(s) + O_2(g) \rightarrow SiO_2(s)$

Just as with metals and metalloids, the product of the reaction of oxygen with a **nonmetal** is an oxide. With nonmetals, the oxides formed are sometimes gaseous. For example, elemental sulfur (**S**₈) reacts with oxygen to form sulfur dioxide, a poisonous gas that can be used as a food additive to sterilize dried fruit and wine.

$\mathbf{S}_{8}(s) + \mathbf{8O}_{2}(g) \rightarrow \mathbf{8SO}_{2}(g)$

Whether a metal oxide compound decomposes easily depends largely on the *strength of the metal-oxygen bonds*. If the bonds are *strong*, a *large amount of energy* (high temperature) is required to break them. If the bonds are *relatively weak*, a *smaller amount of energy* (lower temperature) is required to break them, i.e. the temperature at which the compound decomposes is relatively low.

Procedure

PART 1: DECOMPOSITION OF OXYGEN-CONTAINING COMPOUNDS

Your instructor will demonstrate this part.

Each oxide listed below will be heated using a Bunsen burner. As your instructor presents each case, observe and describe the physical properties of these oxides prior to heating. Physical properties include *color, texture,* and *state of matter*. Thoroughly observe any changes that occur. You must record whether or not there is any change, along with any of the changes in color, texture, gas formation, melting, etc., all in the data sheet below.

Name of Oxide	Chemical Formula	Physical Properties	Observations upon heating
Sodium Nitrate	NaNO₃		
Silicon Dioxide	SiO ₂		
Lead(IV) Oxide	PbO2		
Magnesium Oxide	MgO		
Potassium Chlorate	KClO ₃		

In the space below, write the chemical reaction equations <u>for each substance that</u> <u>released oxygen</u>:

PART 2: COLLECTION OF OXYGEN OVER WATER

You will be collecting **3** bottles of oxygen which will be produced by the *decomposition of hydrogen peroxide* (H_2O_2). The following equation expresses the chemical change that occurs:

 $2 H_2O_2 (I) \xrightarrow{} O_2 (g) + 2H_2O (I)$ MnO₂ (catalyst)

The oxygen is collected by a method known as *downward displacement of water*. The gas is conducted from a generator to an inverted bottle of water through a water-filled trough. The oxygen, which is only very slightly soluble in water, rises and forces the water out.

The rate of decomposition is accelerated using a catalyst, Manganese Dioxide.

PROCEDURE:

Your apparatus will consist of a collecting trough that connects via rubber tubing to a generator bottle with a thistle tube. The reaction will take place in the generator bottle and then the oxygen gas will move through the tubing into the collecting trough where the collection bottles will be.

- 1. Partly fill the collecting trough with tap water, making sure the shelf is *under* the water level.
- Fill 3 collection bottles with tap water. Use the accompanying glass plate to cover a bottle, invert the bottle and place it into the trough upside down. Repeat with the remaining bottles. * * *Before moving on to the next step*, make sure each bottle is *entirely* filled with water upside down in the trough.
- 3. The *generator bottle* already contains the catalyst (**MnO**₂). Once your instructor confirms your setup, *AND* that you are wearing lab gloves, you will receive 40 mL of hydrogen peroxide (**H**₂**O**₂). When instructed to begin, add your H₂O₂ **s I o w I y**.
- 4. Once bubbles start forming in the trough, *wait a few seconds* for all the air to be purged from the generator bottle and tubing.
- 5. Slide one of the collection bottles over the hole that the bubbles are passing through and hold the bottle down with one hand. As oxygen enters the bottle, it will displace the water, and the bottle will become buoyant. This could cause it to flip over, and you will have to start over with that bottle if so.
- 6. Once the bottle fills *almost* entirely with oxygen, you will use one of the bottle lids to seal all of the gas in the bottle. You want to have a small amount of water (2-3 mL) in the bottle before you cap it. The way to do this is lift the bottle straight up without turning it, slide the lid underneath, and close it. *Only after sealing the bottle* will you set it aside. Repeat **5** and **6** with the remaining two bottles.
 - * * * WAIT for your instructor to indicate when to start Part 3, which requires all **3** bottles.

PART 3: PREPARATION OF OXIDES, FORMATION OF ACIDS AND BASES FROM OXIDES

As stated earlier, oxygen is a *reactive element*. To demonstrate this idea we will use the oxygen collected in Part 2 to react with elemental *sulfur*, *phosphorus*, and *magnesium*. In each case, the reaction will produce an <u>oxide</u>.

PROCEDURE

- 1. Label one collection bottle "**sulfur**." Label the second bottle "**phosphorus**." Your instructor will use the third bottle for class demonstrations.
- 2. Bring the bottle labeled "sulfur" to a fume hood. Your instructor will place a small amount of sulfur in a deflagrating spoon. Using the Bunsen burner, the sulfur will be heated until it starts to burn. <u>Quickly</u> open the bottle containing oxygen, and <u>quickly</u> place the spoon over the bottle opening. Allow the sulfur to burn in the oxygen-enriched atmosphere. Replace the lid onto the bottle. Shake the solution in the bottle. Take the bottle back to your workstation.
- 3. Repeat the same procedure (step 2) using phosphorus and the bottle labeled "phosphorus."
- 4. Your instructor will place **two drops** of *universal indicator solution* (UIS) in some of each "sulfur" and "phosphorus" bottle. Record the resulting color (ROYGBIV) of the solution in each bottle.
- 5. Your instructor will demonstrate the remaining three elements.

Record Part 3 Observations Below:

Oxide	Color of Solution with <u>UIS</u>	Color of solution with Phenolphthalein
SO ₂		
P ₂ O ₅		
MgO		
CO2		
Fe ₂ O ₃		

In the space below, write balanced equations for the reaction of each element with OXYGEN.

Now, <u>write balanced equations</u> for the reaction of each <u>oxide</u> with <u>WATER</u>:

Based on your observations, suggest a generalization for acids and bases:

- a) Oxides of *metals*, if dissolved in water, produce ______
- b) Oxides of *nonmetals*, if dissolved in water, produce ______